

## Product datasheet for **GA107221**

### Human IFI44 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	IFI44 CRISPRa kit - CRISPR gene activation of human interferon induced protein 44
Symbol:	IFI44
UniProt ID:	<u>Q8TCB0</u>
Synonyms:	MTAP44; p44; TLDC5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)