

## Product datasheet for **GA107225**

### Human BIG1 (ARFGEF1) activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | ARFGEF1 CRISPRa kit - CRISPR gene activation of human ADP ribosylation factor guanine nucleotide exchange factor 1                                                                                                                                                                                      |
| Symbol:       | ARFGEF1                                                                                                                                                                                                                                                                                                 |
| Components:   | GA107225G1, BIG1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA107225G2, BIG1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA107225G3, BIG1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q9Y6D6</a>                                                                                                                                                                                                                                                                                  |
| Synonyms:     | ARFGEF1; BIG1; P200                                                                                                                                                                                                                                                                                     |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)