

## Product datasheet for **GA107271**

### Human TGN46 (TGOLN2) activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | TGOLN2 CRISPRa kit - CRISPR gene activation of human trans-golgi network protein 2 |
| Symbol:       | TGOLN2   |
| UniProt ID:   | <a href="#">O43493</a>   |
| Synonyms:     | TGN38; TGN46; TGN48; TGN51; TTGN2  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)            |

#### Product images:



[View online >](#)