

Product datasheet for **GA107281**

Human PDPN activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | PDPN CRISPRa kit - CRISPR gene activation of human podoplanin |
| Symbol: | PDPN |
| UniProt ID: | Q86YL7 |
| Synonyms: | AGGRUS; GP36; Gp38; GP40; HT1A-1; OTS8; PA2.26; T1A; T1A-2; T1A2; T11A |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:

