

Product datasheet for **GA107318**

Human **DCTN6** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	DCTN6 CRISPRa kit - CRISPR gene activation of human dynactin subunit 6
Symbol:	DCTN6
Components:	GA107318G1, DCTN6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107318G2, DCTN6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107318G3, DCTN6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O00399</u>
Synonyms:	p27; WS-3; WS3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)