

Product datasheet for **GA107335**

Human TNRC5 (CNPY3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CNPY3 CRISPRa kit - CRISPR gene activation of human canopy FGF signaling regulator 3
Symbol:	CNPY3
UniProt ID:	<u>Q9BT09</u>
Synonyms:	CAG4A; ERDA5; PRAT4A; TNRC5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)