

Product datasheet for **GA107341**

Human TBR1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TBR1 CRISPRa kit – CRISPR gene activation of human T-box brain transcription factor 1
Symbol:	TBR1
Components:	GA107341G1, TBR1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107341G2, TBR1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107341G3, TBR1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q16650</u>
Synonyms:	TBR-1; TES-56
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

