

Product datasheet for **GA107350**

Human p23 (PTGES3) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | PTGES3 CRISPRa kit - CRISPR gene activation of human prostaglandin E synthase 3 |
| Symbol: | PTGES3 |
| UniProt ID: | Q15185 |
| Synonyms: | cPGES; P23; TEBP |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)