

Product datasheet for **GA107415**

Human SPINT3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SPINT3 CRISPRa kit - CRISPR gene activation of human serine peptidase inhibitor, Kunitz type 3
Symbol:	SPINT3
UniProt ID:	P49223
Synonyms:	HKIB9
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)