

Product datasheet for **GA107441**

Human Tetraspanin 9 (TSPAN9) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSPAN9 CRISPRa kit - CRISPR gene activation of human tetraspanin 9
Symbol:	TSPAN9
Components:	GA107441G1, Tetraspanin 9 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107441G2, Tetraspanin 9 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107441G3, Tetraspanin 9 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O75954</u>
Synonyms:	NET-5; NET5; PP1057
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)