

Product datasheet for **GA107445**

Human CD300C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CD300C CRISPRa kit - CRISPR gene activation of human CD300c molecule
Symbol:	CD300C
UniProt ID:	Q08708
Synonyms:	CLM-6; CMRF-35; CMRF-35A; CMRF35; CMRF35-A1; CMRF35A; CMRF35A1; IGSF16; LIR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)