

Product datasheet for **GA107452**

Human SMR3B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SMR3B CRISPRa kit - CRISPR gene activation of human submaxillary gland androgen regulated protein 3B
Symbol:	SMR3B
UniProt ID:	P02814
Synonyms:	P-B; PBII; PRL3; PROL3; SMR1B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)