

Product datasheet for **GA107467**

Human YIF1A activation kit by CRISPRa

Product data:

Product Type: CRISPRa Kits

Description: YIF1A CRISPRa kit - CRISPR gene activation of human Yip1 interacting factor homolog A, membrane trafficking protein

Symbol: YIF1A


Components: GA107467G1, YIF1A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)
GA107467G2, YIF1A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)
GA107467G3, YIF1A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)
1 CRISPRa-Enhancer vector, SKU GE100056
1 CRISPRa scramble vector, SKU GE100077

UniProt ID: [Q95070](#)

Synonyms: 54TM; FinGER7; YIF1; YIF1P

Format: 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)