


Product datasheet for **GA107477**

Human TXNL4A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TXNL4A CRISPRa kit - CRISPR gene activation of human thioredoxin like 4A
Symbol:	TXNL4A
Components:	GA107477G1, TXNL4A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107477G2, TXNL4A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107477G3, TXNL4A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	P83876
Synonyms:	BMKS; DIB1; DIM1; SNRNP15; TXNL4; U5-15kD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)