

Product datasheet for **GA107491**

Human PC4 (SUB1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SUB1 CRISPRa kit - CRISPR gene activation of human SUB1 regulator of transcription
Symbol:	SUB1
UniProt ID:	P53999
Synonyms:	p14; P15; PC4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)