

Product datasheet for **GA107492**

Human SMPDL3A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SMPDL3A CRISPRa kit - CRISPR gene activation of human sphingomyelin phosphodiesterase acid like 3A
Symbol:	SMPDL3A
UniProt ID:	Q92484
Synonyms:	ASM3A; ASML3a; yR36GH4.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)