


## Product datasheet for **GA107507**

### Human SMAP (C11orf58) activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                         |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| Description:  | C11orf58 CRISPRa kit - CRISPR gene activation of human chromosome 11 open reading frame 58                                                                                                                                                                                                              |
| Symbol:       | C11orf58                                                                                                                                                                                                                                                                                                |
| Components:   | GA107507G1, SMAP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA107507G2, SMAP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA107507G3, SMAP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q00193</a></u>                                                                                                                                                                                                                                                                           |
| Synonyms:     | chromosome 11 open reading frame 58; DKFZp686N0232; IMAGE145052; MGC117265; small acidic protein; SMAP                                                                                                                                                                                                  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |

#### Product images:



[View online »](#)