

Product datasheet for **GA107508**

Human KDEL Receptor (KDEL1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KDEL1 CRISPRa kit - CRISPR gene activation of human KDEL endoplasmic reticulum protein retention receptor 1
Symbol:	KDEL1
UniProt ID:	P24390
Synonyms:	ERD2; ERD2.1; HDEL; PM23
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)