

Product datasheet for **GA107509**

Human **SF3A3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SF3A3 CRISPRa kit - CRISPR gene activation of human splicing factor 3a subunit 3
Symbol:	SF3A3
Components:	GA107509G1, SF3A3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107509G2, SF3A3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107509G3, SF3A3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q12874</u>
Synonyms:	PRP9; PRPF9; SAP61; SF3a60
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)