

Product datasheet for **GA107518**

Human **SERINC3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SERINC3 CRISPRa kit - CRISPR gene activation of human serine incorporator 3
Symbol:	SERINC3
Components:	GA107518G1, SERINC3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107518G2, SERINC3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107518G3, SERINC3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q13530
Synonyms:	AIGP1; DIFF33; SBB199; TDE; TDE1; TMS-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)