


Product datasheet for **GA107520**

Human PNRC1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PNRC1 CRISPRa kit - CRISPR gene activation of human proline rich nuclear receptor coactivator 1
Symbol:	PNRC1
Components:	GA107520G1, PNRC1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107520G2, PNRC1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107520G3, PNRC1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q12796</u>
Synonyms:	B4-2; PNAS-145; PROL2; PRR2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)