

Product datasheet for GA107525

Human STIP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	STIP1 CRISPRa kit - CRISPR gene activation of human stress induced phosphoprotein 1
Symbol:	STIP1
UniProt ID:	P31948
Synonyms:	HEL-S-94n; HOP; IEF-SSP-3521; P60; STI1; STI1L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)