

Product datasheet for **GA107526**

Human IFI44L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IFI44L CRISPRa kit - CRISPR gene activation of human interferon induced protein 44 like
Symbol:	IFI44L
UniProt ID:	<u>Q53G44</u>
Synonyms:	C1orf29; GS3686
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)