

Product datasheet for **GA107546**

Human LILRB5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LILRB5 CRISPRa kit - CRISPR gene activation of human leukocyte immunoglobulin like receptor B5
Symbol:	LILRB5
Components:	GA107546G1, LILRB5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107546G2, LILRB5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107546G3, LILRB5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q75023</u>
Synonyms:	CD85C; LIR-8; LIR8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)