


Product datasheet for **GA107547**

Human **SLC38A3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC38A3 CRISPRa kit - CRISPR gene activation of human solute carrier family 38 member 3
Symbol:	SLC38A3
UniProt ID:	Q99624
Synonyms:	G17; NAT1; SN1; SNAT3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)