

Product datasheet for **GA107557**

Human LILRB4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LILRB4 CRISPRa kit - CRISPR gene activation of human leukocyte immunoglobulin like receptor B4
Symbol:	LILRB4
UniProt ID:	Q8NHJ6
Synonyms:	CD85K; HM18; ILT3; LILRB5; LIR-5; LIR5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)