

Product datasheet for **GA107565**

Human KDELR3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KDELR3 CRISPRa kit - CRISPR gene activation of human KDEL endoplasmic reticulum protein retention receptor 3
Symbol:	KDELR3
UniProt ID:	O43731
Synonyms:	ERD2L3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)