


Product datasheet for **GA107574**

Human LILRA1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LILRA1 CRISPRa kit - CRISPR gene activation of human leukocyte immunoglobulin like receptor A1
Symbol:	LILRA1
Components:	GA107574G1, LILRA1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107574G2, LILRA1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107574G3, LILRA1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q75019</u>
Synonyms:	CD85I; LIR-6; LIR6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)