

Product datasheet for **GA107575**

Human LILRB3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LILRB3 CRISPRa kit - CRISPR gene activation of human leukocyte immunoglobulin like receptor B3
Symbol:	LILRB3
UniProt ID:	O75022
Synonyms:	CD85A; HL9; ILT-5; ILT5; LILRA6; LIR-3; LIR3; PIRB
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)