

## Product datasheet for **GA107576**

### Human LILRA3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	LILRA3 CRISPRa kit - CRISPR gene activation of human leukocyte immunoglobulin like receptor A3
Symbol:	LILRA3
UniProt ID:	<u><a href="#">Q8N6C8</a></u>
Synonyms:	CD85E; HM31; HM43; ILT-6; ILT6; LIR-4; LIR4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)