

## Product datasheet for **GA107582**

### Human RIP3 (RIPK3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	RIPK3 CRISPRa kit - CRISPR gene activation of human receptor interacting serine/threonine kinase 3
Symbol:	RIPK3
UniProt ID:	<u><a href="#">Q9Y572</a></u>
Synonyms:	RIP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)