

## Product datasheet for **GA107588**

### Human TRF41 (PAPD7) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TENT4A CRISPRa kit - CRISPR gene activation of human terminal nucleotidyltransferase 4A
Symbol:	TENT4A
UniProt ID:	<a href="#">Q5XG87</a>
Synonyms:	LAK-1; LAK1; POLK; POLS; TRF4; TRF4-1; TRF41; TUTASE5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)