

Product datasheet for **GA107590**

Human SLC35D2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC35D2 CRISPRa kit - CRISPR gene activation of human solute carrier family 35 member D2
Symbol:	SLC35D2
Components:	GA107590G1, SLC35D2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107590G2, SLC35D2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107590G3, SLC35D2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q76EJ3
Synonyms:	hfr; HFRC1; SQV7L; UGTrel8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)