


Product datasheet for **GA107601**

Human DUS4L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DUS4L CRISPRa kit - CRISPR gene activation of human dihydrouridine synthase 4 like
Symbol:	DUS4L
Components:	GA107601G1, DUS4L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107601G2, DUS4L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107601G3, DUS4L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O95620</u>
Synonyms:	DUS4; PP35
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)