

Product datasheet for **GA107603**

Human CNTRL activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CNTRL CRISPRa kit - CRISPR gene activation of human centriolin
Symbol:	CNTRL
Components:	GA107603G1, CNTRL gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107603G2, CNTRL gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107603G3, CNTRL gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q7Z7A1</u>
Synonyms:	bA165P4.1; CEP1; CEP110; FAN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)