

Product datasheet for **GA107609**

Human PL6 (TMEM115) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM115 CRISPRa kit – CRISPR gene activation of human transmembrane protein 115
Symbol:	TMEM115
Components:	GA107609G1, PL6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107609G2, PL6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107609G3, PL6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q12893
Synonyms:	PL6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

