

Product datasheet for **GA107611**

Human TOPBP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TOPBP1 CRISPRa kit - CRISPR gene activation of human DNA topoisomerase II binding protein 1
Symbol:	TOPBP1
UniProt ID:	Q92547
Synonyms:	TOP2BP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)