

Product datasheet for **GA107661**

Human PWP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PWP1 CRISPRa kit - CRISPR gene activation of human PWP1 homolog, endonuclein
Symbol:	PWP1
UniProt ID:	<u>Q13610</u>
Synonyms:	IEF-SSP-9502
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)