

## Product datasheet for **GA107676**

### Human AP4S1 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	AP4S1 CRISPRa kit - CRISPR gene activation of human adaptor related protein complex 4 subunit sigma 1
Symbol:	AP4S1
UniProt ID:	<u><a href="#">Q9Y587</a></u>
Synonyms:	AP47B; CLA20; CLAPS4; CPSQ6; SPG52
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)