

## Product datasheet for **GA107690**

### Human PSIP1 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PSIP1 CRISPRa kit - CRISPR gene activation of human PC4 and SFRS1 interacting protein 1
Symbol:	PSIP1
UniProt ID:	<a href="#">O75475</a>
Synonyms:	DFS70; LEDGF; p52; p75; PAIP; PSIP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)