

Product datasheet for **GA107693**

Human Unrip (STRAP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	STRAP CRISPRa kit - CRISPR gene activation of human serine/threonine kinase receptor associated protein
Symbol:	STRAP
UniProt ID:	Q9Y3F4
Synonyms:	MAWD; PT-WD; UNRIP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)