

## Product datasheet for **GA107714**

### Human SEC23IP activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SEC23IP CRISPRa kit - CRISPR gene activation of human SEC23 interacting protein
Symbol:	SEC23IP
UniProt ID:	<u><a href="#">Q9Y6Y8</a></u>
Synonyms:	MSTP053; P125; P125A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)