


Product datasheet for **GA107716**

Human SUPT16H activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SUPT16H CRISPRa kit - CRISPR gene activation of human SPT16 homolog, facilitates chromatin remodeling subunit
Symbol:	SUPT16H
UniProt ID:	Q9Y5B9
Synonyms:	CDC68; FACTP140; SPT16; SPT16/CDC68
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)