

Product datasheet for **GA107719**

Human DNA Polymerase iota (POLI) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	POLI CRISPRa kit - CRISPR gene activation of human DNA polymerase iota
Symbol:	POLI
UniProt ID:	Q9UNA4
Synonyms:	RAD3OB; RAD30B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)