

Product datasheet for **GA107727**

Human Gemin 3 (DDX20) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DDX20 CRISPRa kit - CRISPR gene activation of human DEAD-box helicase 20
Symbol:	DDX20
Components:	GA107727G1, Gemin 3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107727G2, Gemin 3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107727G3, Gemin 3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9UHI6
Synonyms:	DP103; GEMIN3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)