

Product datasheet for **GA107752**

Human GPR44 (PTGDR2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PTGDR2 CRISPRa kit - CRISPR gene activation of human prostaglandin D2 receptor 2
Symbol:	PTGDR2
UniProt ID:	<u>Q9Y5Y4</u>
Synonyms:	CD294; CRTH2; DL1R; DP2; GPR44
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)