

Product datasheet for **GA107762**

Human PXMP4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PXMP4 CRISPRa kit - CRISPR gene activation of human peroxisomal membrane protein 4
Symbol:	PXMP4
Components:	GA107762G1, PXMP4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107762G2, PXMP4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107762G3, PXMP4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y6I8</u>
Synonyms:	PMP24
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)