

Product datasheet for **GA107779**

Human PNK (PNKP) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | PNKP CRISPRa kit - CRISPR gene activation of human polynucleotide kinase 3'-phosphatase |
| Symbol: | PNKP |
| UniProt ID: | Q96T60 |
| Synonyms: | AOA4; EIEE10; MCSZ; PNK |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)