

Product datasheet for **GA107784**

Human CD300 antigen (CD300A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CD300A CRISPRa kit - CRISPR gene activation of human CD300a molecule
Symbol:	CD300A
UniProt ID:	<u>Q9UGN4</u>
Synonyms:	CLM-8; CMRF-35-H9; CMRF-35H; CMRF35-H; CMRF35-H9; CMRF35H; CMRF35H9; IGSF12; IRC1; IRC1/IRC2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)