


Product datasheet for **GA107811**

Human TWF2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TWF2 CRISPRa kit - CRISPR gene activation of human twinfilin actin binding protein 2
Symbol:	TWF2
Components:	GA107811G1, TWF2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107811G2, TWF2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107811G3, TWF2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6IBS0</u>
Synonyms:	A6r; A6RP; MSTP011; PTK9L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)